

This is a feature list to compare our solution with competitive products.

	BIOBSERVE	Competitor	Competitor
Recording technique			
No software required on test computer/device	✓		
Losless monitor recording with the full resolution	✓		
All operating systems can be recorded	✓		
Any VGA signal can be recorded	✓		
No picture-in-picture	✓		
Multiple video and audio recordings possible	✓		
Multiple monitor recordings possible	✓		
Easy to set up	✓		
Spectator software			
Configuration			
Participant database (freely definable)	✓		
Projet database (up to five levels)	✓		
Unlimited behaviors	✓		
Unlimited modifiers for each behavior	✓		
Multi key codes	✓		
Task lists	✓		
Task rating (definable)	✓		
Session composer (combine task and behavior lists)	✓		
Session calendar	✓		
Easy use of different set up configurations	✓		
Backup and restore tool	✓		
Acquisition			
Coding with keys or mouse	✓		
General comments	✓		
Comments attached to behaviors	✓		
Real time data visualization (table and time chart)	✓		
Undo command	✓		
New behaviors configurable during the session	✓		
Multi Observer plug in available	✓		
Play & Edit			
Video and data replay	✓		
All acquired data editable	✓		
Frame by frame analysis	✓		
Fast forward and slow motion	✓		
Analysis			
Highlighted Video	✓		
Presentation tool (show data and videos)	✓		
Report generator	✓		
Pattern analysis	✓		
Statistical Analysis of Tasks and Behaviors	✓		
Reliability check of multiple observers	✓		
Participant filter	✓		
Summaries of multiple test sessions	✓		

If you miss a feature or if a competitor offers a feature that is not listed, please contact us.